

188

Injection of 2 rabbits with Linear 78-109 peptide
 no 86-87 from TGF β .

2 " with cyclized 78-109 "
 no 88-89

Procedure: 2 mg / rabbit of each peptide

- 2 mg peptide were dissolved in 250 μ l DAW
- 0.5 mg Methylated BSA were ~~added to the solution to couple the peptide~~
~~added to the solution to couple the peptide~~
 (Sigma A1009: MeBSA)
- vortex 30 "
- Added 250 μ l Freund's adjuvant complete
- mixed 1 h with homogenizer
- 0.5 ml solution was injected in each rabbit

Note: Both peptides were difficult to dissolve.
 " " were not HPLC purified

Injection of 2 rabbits: Linear TGF- β peptide

" " Cyclic "

" " PG peptide

Procedure

except Freund's incomplete adjuvant was used.

For Wayne: KLM-peptide was already supplied
I only mixed to F.I. 0.8% mix 1h 4°C
inject 0.5 ml / rabbit

Rabbit MF 1285 for PG AUFI 1385
1285

ANIMAL PROCEDURE REQUEST

(For procedures to be performed by Animal Facility personnel)
Return this form to the Animal Facility Office.

LJCRF

TODAY'S DATE: _____

| | | | | | |
|-------------------------|-------------------|-------|-----|---------|--|
| Principal Investigator | E. RUSLAHTI | Phone | 232 | Lab No. | |
| Investigator/Technician | LUCIA R. LANGVINO | Phone | | Lab No. | |

This procedure relates to the project outlined on Animal Usage Form No. _____

Date or Dates on which this procedure is required (please allow 7 days for scheduling):

| | | | |
|---|-------------|--|---|
| 1 | 50 ml | (If more than 8 procedures are required, attach a second form) | 5 |
| 2 | " | | 6 |
| 3 | Euth 100 ml | | 7 |
| 4 | | | 8 |

Describe the procedure(s) required:

Vol. of blood to be collected each procedure: _____ ml.

8 Injection 6 rabbits 2 PC
2 Tiff linear probe
2 cyclic u